NOTA TÉCNICA CONJUNTA COVISA/CAB N° 01/2024

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English

About the test

[Introduction] Dengue viruses, transmitted by *Ae. aegypti* and *A. albopictus* mosquitoes, are widely distributed throughout the tropical and subtropical areas of the world. Dengue is an important arthropod-borne viral disease due to its high morbidity and mortality in humans. There are four known distinct serotypes of dengue virus, 2, 3 and 4. In children, infection is often self-limited or causes a self-limited febrile disease. However, if the patient is infected with a different serotype, risk of a more severe disease increases, such as dengue hemorrhagic fever or dengue shock syndrome. NS1 antigen is a highly-conserved glycoprotein that is present at high concentrations in the sera of dengue secondary infections, during the early clinical phase of the disease. NS1 antigen is found in specimens from primary or secondary dengue infected patients between 1 and 7 days after onset of fever. IgM usually becomes detectable only 10 to 14 days after the onset of illness in cases of primary dengue infection and only 4 to 5 days after onset of illness in secondary infections. In primary infections, IgG appears after the 14th day and persists for 18 months. In secondary infections, IgG levels rise to 1 - 2 days after the onset of symptoms and an IgM response is observed 20 days after infection.

[Principle] The SB BLOLINE Dengue Duo Rapid Test contains two in vitro immunochromatography assays for analysis of dengue infection in human serum, plasma or whole blood (see table; Dengue NS1 Ag assay, right side; Dengue IgG/IgM assay). The Dengue NS1 antigen Rapid Test is a highly sensitive and specific one step assay for the qualitative detection of dengue virus NS1 antigen. When added to the specimen well, NS1 antigens present in the specimen reacts with mouse monoclonal anti-NS1 Ag-colditoid gold conjugates, forming a complex of antibodies and antigens. As this complex migrates along the test device, by chromatography, it is captured by immobilized antibodies conjugated to latex particles (anti-NS1 antigen Test Line). A second line is used for procedural control and always appear adjacent to C (Control) line if the test procedure is performed correctly and the active ingredients of main components on the strip are working. When dengue antibodies are present in the specimen, both the Control (C) and Test (T) Lines appear in the result window, indicating a positive result. Only the Control Line (C) appears for specimens that do not contain NS1 antigen or contain NS1 antigen below the detectable level. The Dengue IgG/IgM Rapid Test is a qualitative assay for the simultaneous detection of IgG and IgM antibodies to dengue virus. This test detects antibodies for all four dengue serotypes by using a mixture of recombinant dengue envelope proteins. When added to the specimen well, anti-dengue IgG and IgM present in the specimen react with the recombinant dengue virus envelope proteins-colditoid gold conjugates, forming a complex of antibodies and antigens. This complex migrates along the test device, by chromatography, to a capture strip with two colored Test Lines adjacent to "G" (Dengue IgG Test Line) and/or "M" (Dengue IgM Test Line). A third line is used for procedural control and always appear adjacent to C (Control) line if the test procedure is performed correctly and the active ingredients of main components on the strip are working. When dengue antibodies are present in the specimen, both the Control (C) and Test (G and/or M) Lines appear in the result window, indicating a positive result. Only the Control Line (C) appears for specimens that do not contain dengue virus or dengue virus below the detectable level. The Dengue IgG/IgM Rapid Test is a qualitative assay for the differential detection of IgG and IgM antibodies to dengue virus. This test detects antibodies for all four dengue serotypes by using a mixture of recombinant dengue envelope proteins. When added to the specimen well, anti-dengue IgG and IgM present in the specimen react with the recombinant dengue virus envelope proteins-colditoid gold conjugates, forming a complex of antibodies and antigens. This complex migrates along the test device, by chromatography, to a capture strip with two colored Test Lines adjacent to "G" (Dengue IgG Test Line) and/or "M" (Dengue IgM Test Line). A third line is used for procedural control and always appear adjacent to C (Control) line if the test procedure is performed correctly and the active ingredients of main components on the strip are working. When dengue antibodies are present in the specimen, both the Control (C) and Test (G and/or M) Lines appear in the result window, indicating a positive result. Only the Control Line (C) appears for specimens that do not contain dengue virus or dengue virus below the detectable level.

[Intended use] The SB BLOLINE Dengue Duo Rapid Test is an in vitro immunochromatography, one step assay designed to detect dengue virus NS1 antigen, IgG and IgM antibodies to dengue virus in human serum, plasma or whole blood. This test is intended for professional use, to aid in the presumptive diagnosis of primary and secondary dengue infections. The SB BLOLINE Dengue Duo Rapid Test is not a confirmatory test. Specific alternative diagnostic methods must be used to obtain a confirmation of dengue virus infection. This may include isolation of virus, antigen detection in fixed tissues, PCR and serological testing such as haemagglutination inhibition.

Materials provided and active ingredients of main components

1. The test kit contains the following items:

Component	Cat. No. 11K145	Cat. No. 11K146
Dengue NS1 Ag and Dengue IgG/IgM combo cassette, individually packaged including desiccants sachet	10	25
Assay diluent for Dengue IgG/IgM test	1.3ml	1.5ml
Capillary pipettes for dengue IgG/IgM test (10 µl)	10	25
Disposable droppers for dengue NS1 Ag test	10	25
Instructions for use	1	1

2. Active ingredients of main components

1. **Dengue NS1 Ag:** Gold conjugate : Mouse monoclonal anti-dengue NS1 – gold colloid (0.2,0.3,0.5 µg). Test Line G : Mouse monoclonal anti-dengue NS1 (0.2,0.6,1 µg). Control Line : Goat anti-mouse IgG (0.2,0.3,0.5 µg). Control Line : Rabbit anti-mouse IgG (3.5,2.5 µg).

2. **Dengue IgG/IgM:**

1. **1st strip:** gold conjugate : Recombinant dengue virus envelope protein – gold colloid (1.0,1.2 µg). Test Line G : Mouse monoclonal anti-human IgG (5=1 µg). Test Line M: Mouse monoclonal anti-human IgM (5=1 µg). Control Line : Rabbit anti-mouse IgG (3.5,2.5 µg).

2. **Assay diluent:** 1.3ml, Sodium azide (0.01 w/v)

Kit storage and stability

- Store the test kit in a cool, dry and humidify and must be used and stored properly.
- Store the test kit at between 1 °C and 30 °C. Do not freeze the kit or its components or expose them to high temperatures.
- Before using the test, check the humidity indicator on the desiccants on the pouch for color change. Discard the device if the humidity indicator on the device is yellow.
- Perform the test immediately after removing the test device from the foil pouch, to avoid exposure to moisture.

3. Do not use the kit or its components beyond the expiration date.

4. Do not use the test device if the pouch is damaged or the seal is broken.

Warnings

- For a viral diagnostic use only. Do not reuse test device and kit components.
- The instructions must be followed exactly to achieve accurate results. Any individual performing an assay using this product must be trained in its use and must be experienced in laboratory procedures.
- Do not eat or smoke while handling specimens.
- Wear protective gloves while handling specimens and wash hands thoroughly after procedures.
- Avoid splashing or aerosol formation of specimens and assay diluent.
- Clean up spills thoroughly using appropriate disinfectant.
- Decomtaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. specimen applicator, test device) in a biohazard container as if they were infectious waste.

Specimen collection and storage

Whole blood

- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).

Caution: Do not use capillary blood.

- Hold specimens and test immediately if tested, they must be refrigerated at 2 - 8 °C.
 - If stored at 2 - 8 °C, the blood specimens must be tested within 3 days.
 - Do not use a blood specimen stored for more than 3 days; it can cause a non-specific reaction.
 - Bring blood specimens to room temperature (15 - 30 °C) prior to use.
- Plasma or serum**
- Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge the tube to get a plasma specimen.
 - [Serum]** Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulants), allow the blood to coagulate for 30 minutes and then centrifuge the tube to get a serum specimen.
 - If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C for storage period longer than 2 weeks, freezing below -15 °C is required. Specimens must be brought to room temperature (15 °C - 30 °C) prior to use.
 - Plasma or serum specimens containing a precipitate may yield inconsistent test results. Use standard laboratory techniques to clarify such specimens prior to testing.

Precautions

- Repaired gloves should only be used if the specimens should be avoided.
- Anticoagulants including heparin, EDTA, and citrate do not affect the test result. Use of other anticoagulants may not have been validated. They are not affect the test result.
- Use a new disposable dropper or pipette for test each specimen in order to avoid cross-contamination of specimens, which could produce erroneous results.

Test procedure (Refer to figure)

- Allow all kit components and specimen to reach room temperature (Between 15 ° and 30 °C) prior to opening and testing.
- Remove the test device from the foil pouch, and place it on a flat dry surface. Label the device with a patient identifier. Perform the test immediately to avoid exposing test device to moisture.

Dengue NS1 Ag					
1. With a disposable dropper provided, dispense 3 drops (about 100 µl) of specimen into the specimen well.	Primary	89			
	Intermediate	2			
	Secondary	1			
2. With 10 µl capillary pipette provided, draw specimen to black line (10 µl) then dispense into the square specimen well.	Primary	81	245 ¹	206	39
	Intermediate	1			
	Secondary	1			
3. Dispense 3 drops (about 90 - 120 µl) of assay diluent into the round assay diluent well.	Primary	8			
	Intermediate	1			
	Secondary	1			
4. Interpret test results at 15 - 20 minutes.	Primary	8			
	Intermediate	1			
	Secondary	1			
Other febrile illness	Primary	47	1	46	
	Secondary	1			

¹150 panels were collected less than 3 days of illness and 89 panels were collected more than 3 days of illness.

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3. **Study Site 3**
Total 259 patients serum specimens collected from Sri Lanka were tested on the SD BIOLINE Dengue Duo. The following 4 potential cross-reacting pathogens had no effect on test results of SD BIOLINE Dengue Duo. Even though there are no cross reactivity in the results of internal study, the possibility can't be excluded completely.

Disease	Dengue IgM Positive/Total	Dengue IgG Positive/Total	Dengue NS1 Ag Positive/Total
Japanese Encephalitis	0/25	0/25	0/25
Yellow Fever	0/25	0/25	0/25
Malaria P. falciparum	0/25	0/25	0/25
Malaria P. vivax	0/25	0/25	0/25
Total	0/100	0/100	0/100

7. Reproducibility of SD BIOLINE Dengue Duo test results has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria. The performance of SD BIOLINE Dengue Duo kit was evaluated with potential interfering substances. Hemolytic, lipemic, icteric specimens and those containing rheumatoid factors were investigated. None of the specimens were found to interfere with this test kit.

4. **Study Site 4**
Sensitivity of SD BIOLINE Dengue Duo according to the type and phase of dengue infection
A total of 397 serum specimens collected from Mexico were tested on the SD BIOLINE Dengue Duo in order to evaluate the performance. The specimens were confirmed as having dengue infection by ELISA test. On the basis of NS1 Ag & IgG/IgM test results and clinical date of patients, 310 specimens were diagnosed with dengue virus infection and 87 specimens were negative for all 3 ELISAs. The performance of SD BIOLINE Dengue Duo was evaluated against the overall results of ELISA tests.

Specimen type	Specimen No.	SD BIOLINE Dengue Duo Result (95% CI)	Sensitivity
Primary acute dengue	122	110 12 90.16% (84.47 – 95.86%)	96.63%
	17	11 6 64.71% (39.05 – 90.36%)	
	122	117 5 95.90% (91.97 – 99.83%)	
Secondary acute dengue	49	43 6 87.76% (77.56 – 97.95%)	94.05%
Non dengue	87	9 78 Specificity: 89.66% (82.68 - 96.63%)	

5. **Study Site 5**
A total of 199 serum specimens collected and tested on the SD BIOLINE Dengue Duo at Ministerio de Salud (MINSAL), Argentina. The serum specimens were comprised of 149 positive and 50 negative specimens confirmed by ELISA and RT-PCR. The performance of SD BIOLINE Dengue Duo was evaluated against the overall results of ELISA and RT-PCR.

Specimen Type	Days after symptom onset	No. of specimens	SD BIOLINE Dengue Duo Result (95% CI)	Sensitivity
Positive	<3	49	24 25 90.16% (84.47 – 95.86%)	96.63%
	4-6	50	35 15 95.90% (91.97 – 99.83%)	
	>6	50	50 0 100%	
Negative	50	0	50 0 100%	96.63%
	122	110 12 90.16% (84.47 – 95.86%)		
	17	11 6 64.71% (39.05 – 90.36%)		

6. **Cross-reactivity test with other Flaviviruses mediated and mosquitoes-borne disease**
The following 4 potential cross-reacting pathogens had no effect on test results of SD BIOLINE Dengue Duo. Even though there are no cross reactivity in the results of internal study, the possibility can't be excluded completely.

Disease	Dengue IgM Positive/Total	Dengue IgG Positive/Total	Dengue NS1 Ag Positive/Total
Japanese Encephalitis	0/25	0/25	0/25
Yellow Fever	0/25	0/25	0/25
Malaria P. falciparum	0/25	0/25	0/25
Malaria P. vivax	0/25	0/25	0/25
Total	0/100	0/100	0/100

7. Reproducibility of SD BIOLINE Dengue Duo test results has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria. The performance of SD BIOLINE Dengue Duo kit was evaluated with potential interfering substances. Hemolytic, lipemic, icteric specimens and those containing rheumatoid factors were investigated. None of the specimens were found to interfere with this test kit.

8. **Reproducibility of SD BIOLINE Dengue Duo test results** has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria. The performance of SD BIOLINE Dengue Duo kit was evaluated with potential interfering substances. Hemolytic, lipemic, icteric specimens and those containing rheumatoid factors were investigated. None of the specimens were found to interfere with this test kit.

Product Disclaimer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is and outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

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